



Aqua Creek Products

The Total Solution for Pool and Spa Access

The ABLE™ Aquatic Bathtub Lift Elite

Item# F-25ABLE

The Aquatic Bathtub Lift Elite or ABLE is a wall mounted bathtub lift.

It provides those that need added assistance the ability and security to use their bathtub easily and comfortably.

The manually pivoting seat provides complete control on where you enter the tub and is thin enough to allow you to get as deep into your tub as possible.

Site construction and product assembly is required.



Features:

- 300 lb (136 kg) weight capacity
- Powdercoated stainless steel construction
- Rechargeable battery operation including charger
- Submersible 4 button remote
- Folding seat with seatbelt
- Requires little to no assistance to operate
- Maximum tub height 22"
- Limited lifetime warranty



104 E. CALTON RD. SUITE 110
LAREDO, TX 78041
PH. (956) 723-0770

WWW.SPECIALMOBILITY.COM

