



## Aqua Creek Products

The Total Solution for Pool and Spa Access

# The EZ™ Pool Lift

Item# F-03EZPL

The EZ™ Lift is an affordable, manually operated lift ideal for residential applications. The lift can be used on in ground spas and can also accommodate pool wall heights up to 26 inches (66cm).

The EZ™ Lift is constructed of stainless steel and protected against corrosion with a durable, white epoxy powder coated finish. The lift comes standard with mesh sling and all components to make the lift fully functional.


The EZ™ Lift now comes with a rotation bar to allow attendant greater ease in rotating the unit.

### Features:

- 400 Pound (181 kg) weight capacity
- Manually operated hydraulic ram
- Clears walls up to 26" (66 cm) high
- Durable mesh sling- standard size: 38"x20"
- Adjustable stainless steel chains
- Limited lifetime structural warranty

### Options:

- *Hard seat option - Item #F-0320HSO*
- *Youth sling - Item #F-34EZCS /34"x18" (86x46cm)*
- *Extra long sling - Item #F-037EZLS / 50"x28" (127x71cm)*
- *Standard sling with rigid insert & seat belt - Item #F-036RSS*
- *Additional standard in ground anchor - Item #F-033EZA*
- *Wood deck anchor - Item #F-30600*
- *Surface mount anchor - Item #F-035EZSM*
- *Upgrade Pack - Item #F-03EZPLAP - Includes rigid sling with seat belt, sling bars and cover*
- *Cover - Item #F-440EZC*
- *Hard seat cover - Item #F-440EZHSC*



Shown with optional Surface Mount Anchor



Hard Seat Option



Sling with rigid seat and belt

**PROUDLY SOLD BY:**



104 E. CALTON RD. SUITE 110  
LAREDO, TX 78041  
PH. (956) 723-0770

[WWW.SPECIALMOBILITY.COM](http://WWW.SPECIALMOBILITY.COM)