



Aqua Creek Products
The Total Solution for Pool and Spa Access

The Patriot™ Lift

Item# F-12PPL-HD-AT1

Inspected And Verified to Meet
Applicable Sections of ADA
Guidelines Section 1009.2 By:

wyle
laboratories

A Registered ISO 9001
Laboratory

U.S. Patent No.
5790995



The Patriot Portable Pro Pool™ Lift is the first completely portable aquatic lift that is **independently verified to meet ADA guidelines** and capable of lifting a 450lb* (204kg) load.

This lift requires no anchors and can easily be removed from poolside when not needed. It features a compact, low profile frame that does not intrude into the pool when not being used. A unique sliding mechanism allows the lift to be transported easily on its cart and then lowered onto the pool deck for use without using any fastening hardware.

Features:

- 450 Pound* (204 kg) weight capacity
- Independent 3rd party verified to meet ADA standards
- Rechargeable battery with wall mount charger
- Adjustable footrest & lap belt
- Dual flip-up armrests
- Simple push button operation
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty

Options:

- Upgrade Pack - Item #F-12PPAP
Includes headrest, chest strap, battery & cover
- Headrest - Item #F-422PLH
- Chest strap - Item #F-423CS
- Protective cover - Item #F-12PPC
- Pull out leg rest - Item #F-105LAR
- Cycle Attachment - Item #F-019CA





*450 Pound weight capacity for water drafts up to 6”
400 lb. capacity for other applications.

PROUDLY SOLD BY:



104 E. CALTON RD. SUITE 110
LAREDO, TX 78041
PH. (956) 723-0770

WWW.SPECIALMOBILITY.COM



Push lift to deck side, set brakes, slide into place and it's ready to use.