



Aqua Creek Products

The Total Solution for Pool and Spa Access

The Portable Pro Pool™ Lift

Item# F-004PPPB

The Portable Pro Pool™ Lift is a completely portable aquatic lift that can accommodate loads up to 400lb (181kg). It requires no anchors and can easily be removed from poolside when not needed.

The lift features a compact, low profile frame that does not intrude into the pool when not being used. It's aesthetically pleasing and compact design make it a popular lift.

Features:

- 400 Pound (181 kg) weight capacity
- Rechargeable battery with wall mount charger
- Adjustable footrest & lap belt
- Dual flip-up arm rests
- Simple push button operation
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty

Options:

- *Upgrade Pack- Item #F-004PPPBPB*
Includes: battery, headrest, chest strap and cover
- *Headrest - Item #F-422PLH*
- *Chest strap - Item #F-423CS*
- *Protective cover - Item #F-444PPPC*
- *Cycle Attachment - Item #F-019CA*



PROUDLY SOLD BY:



104 E. CALTON RD. SUITE 110
LAREDO, TX 78041
PH. (956) 723-0770

WWW.SPECIALMOBILITY.COM

U.S. Patent No.
5790995