



Aqua Creek Products

The Total Solution for Pool and Spa Access

The Power EZ™ Pool Lift

Item# F-03PEZ

The Power EZ™ Lift is a great economical choice for individuals needing to gain access to their pool or spa.

The Power EZ™ Lift provides a fully automatic actuator to raise and lower with 360° manual rotation. The lift can accommodate spa walls up to 26 inches (66cm) high and can be set back up to 20 inches (50cm) from the pool edge.

The Power EZ™ Lift now comes with a rotation bar to allow attendant greater ease in rotating the unit.

Features:

- 400 Pound (181 kg) weight capacity
- Rechargeable battery operation with charger
- Clears walls up to 26" (66 cm) high
- Durable mesh sling- standard size: 38"x20"
- Adjustable stainless steel chains
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty

Options:

- *Hard seat option - Item #F-0320HSO*
- *Youth sling - Item #F-34EZCS /34"x18" (86x46cm)*
- *Extra long sling - Item #F-037EZLS / 50"x28" (127x71cm)*
- *Standard sling with rigid insert & seat belt - Item #F-036RSS*
- *Additional standard in ground anchor - Item #F-033EZA*



- *Wood deck anchor - Item #F-30600*
- *Surface mount anchor - Item #F-035EZSM*
- *Upgrade Pack - Item #F-03PEZAP*
Includes battery, rigid sling with seat belt, sling bars and cover
- *Cover - Item #F-440EZC*
- *Hard seat cover - Item #F-440EZHSC*



Sling with rigid seat and belt



Hard Seat Option

PROUDLY SOLD BY:



104 E. CALTON RD. SUITE 110
LAREDO, TX 78041
PH. (956) 723-0770

WWW.SPECIALMOBILITY.COM