



Aqua Creek Products
The Total Solution for Pool and Spa Access

The Pro Pool Lift™

Item# F-004PLB



Inspected And Verified to Meet
Applicable Sections of ADA
Guidelines Section 1009.2 By:

wyle
laboratories

A Nationally Recognized
Test Laboratory

The Pro Pool Lift™ is the first ever aquatic lift to be **verified by an independent third party to be ADA compliant and safe for commercial use.**

The unique design of this lift minimizes the amount of required deck space needed to install and it also does not require the lift to rotate or extend into the water unless being used.

The lift features all stainless steel construction with a protective white epoxy powder coat finish. Ideal for facilities with limited deck space and where aesthetics are of the utmost importance.

Features:

- 450 Pound* (204kg) weight capacity
- Independent 3rd party verified to meet ADA standards
- Unique low profile and small footprint
- Adjustable footrest & lap belt
- Dual flip-up arm rests
- Simple push button operation
- Rechargeable battery with wall mount charger
- Limited lifetime structural warranty

Options:

- *Extended anchor set backs - Item #F-004PLB-23, #F-004PLB-26*
- *Deep water drafts (greater than 9") - Item #F-004PLB-DD12*
- *Reverse configuration (left facing seat) - Item #F-004PLB-R*
- *Upgrade Pack- Item #F-004PLBAP*
Includes battery, headrest, chest strap & cover
- *Headrest - Item #F-422PLH*
- *Chest strap - Item #F-423CS*
- *Pull out leg rest - Item #F-105LAR*
- *Protective cover - Item #F-442PPC*
- *Additional standard anchor - Item #F-433PP4A*
- *Quick attach anchor kit - Item #F-046QAB*
- *Wood deck anchor - Item #F-40220*
- *Transport cart - Item #F-428HT*
- *Portable conversion kit - Item #F-010POR*
- *Cycle attachment - Item #F-019CA*



*450 Pound weight capacity for water drafts up to 6", 400 lb. capacity for other applications.



PROUDLY SOLD BY:



104 E. CALTON RD. SUITE 110
LAREDO, TX 78041
PH. (956) 723-0770

WWW.SPECIALMOBILITY.COM