



Aqua Creek Products

The Total Solution for Pool and Spa Access

The Pro Spa 40 Lift™

Item# F-006PSLB

The Pro Spa 40 Lift™ is a fully powered, ADA compliant* lift designed especially for use with above ground spas and pools with wall heights of 20" up to 40".

The lift provides easy access at the push of a button and is constructed of low carbon stainless steel and protected with a durable white epoxy powder coated finish.

Field measurements are required.

Features:

- 400 Pound (181 kg) weight capacity
- ADA compliant for walls up to 37" high
- Clears walls up to 40" high, 14" thick
- Rechargeable battery operation with charger
- Submersible hand held remote
- Pull out leg rest
- Dual flip up arms
- Limited lifetime structural warranty

Options:

- *Upgrade Pack - Item #F-005SLBAP (must be purchased with lift) Includes extra battery, headrest, chest strap & cover*
- *Headrest - Item #F-422PLH (Can't retrofit)*
- *Chest strap - Item #F-423CS*
- *Wood deck anchor - Item #F-106SLWDP*
- *Protective cover - Item #F-445SLC*



Clears walls up to 40" high!



*ADA Compliant only for walls heights up to 37" high

PROUDLY SOLD BY:



104 E. CALTON RD. SUITE 110
LAREDO, TX 78041
PH. (956) 723-0770

WWW.SPECIALMOBILITY.COM