



Aqua Creek Products

The Total Solution for Pool and Spa Access

The Pro Spa 60 Lift™

Item# F-007EXA

The ADA Compliant* Pro Spa 60 Lift™ is designed to provide access to pools and spas with wall heights of 20 inches and up to 60 inches.

The lift can also be used to accommodate multiple bodies of water that are next to each other and can accommodate patrons up to 300 lbs.

Field measurements are required.

Features:

- 300 Pound (136 kg) weight capacity
- ADA compliant for walls up to 57" high
- Clears walls up to 60" high, 14" thick
- Rechargeable battery operation with charger
- 360 degree fully powered rotation
- Submersible hand held remote
- Pull out leg rest
- Dual flip up arms
- Limited lifetime structural warranty

Options:

- *Upgrade Pack - Item #F-005SLBAP*
Must be purchased with lift - Includes battery, headrest, chest strap & cover
- *Head rest - Item #F-422PLH (Can't retrofit)*
- *Chest strap - Item #F-423CS*
- *Wood deck anchor - Item #F-106SLWDP*
- *Protective cover - Item #F-460PSLC*



Clears walls up to 60" high!



PROUDLY SOLD BY:



104 E. CALTON RD. SUITE 110
LAREDO, TX 78041
PH. (956) 723-0770

WWW.SPECIALMOBILITY.COM

*ADA Compliant only for walls heights up to 57" high