



Aqua Creek Products
The Total Solution for Pool and Spa Access

The Revolution™ Lift

Item# F-700RL



Inspected And Verified to Meet
Applicable Sections of ADA
Guidelines Section 1009.2 By:

wyle
laboratories

A Registered ISO 9001
Laboratory

The Revolution™ Lift is **verified by an independent third party to be ADA compliant.** This lift provides the ability to service multiple bodies of water as well as both above ground spas and in ground pools and spas. Multiple options makes this lift very versatile.

Constructed of all stainless steel materials and protected with a white epoxy powder coated finish, The Revolution™ Lift is the only lift your facility needs.

Features:

- 500 Pound (227 kg) weight capacity
- Rechargeable battery with wall mount charger
- Clears 30" high walls, up to 10" thick
- 360 degree powered rotation in either direction
- Adjustable footrest & lap belt
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty

Options:

- *Upgrade Pack - Item #F-700RLAP*
Includes battery, headrest, chest strap, cover
- *Additional standard anchor - Item #F-710RLA*
- *Headrest - Item #F-722HRR*
- *Chest strap - Item #F-423CS*
- *Protective cover - Item #F-720RLC*
- *Spine board assembly - Item #F-734RSA*
- *Gurney assembly - Item #F-730RSA*
- *Sling seat option - Item #F-706RLSS*
- *Wheelchair attachment - Item #F-705*
- *Transport cart - Item #F-714RTC*
- *Cycle attachment - Item #F-019CA*



PROUDLY SOLD BY:



104 E. CALTON RD. SUITE 110
LAREDO, TX 78041
PH. (956) 723-0770

WWW.SPECIALMOBILITY.COM



Spineboard Option



Gurney Option



Wheelchair Option



Sling Seat Option



Transport Cart