



Aqua Creek Products

The Total Solution for Pool and Spa Access

The Super Power EZ™ Lift

Item# F-190PHBL

The Super Power EZ™ Lift is a modified version of our very popular Power EZ™ Lift to be used on above ground backyard swimming pools.

The Super Power EZ™ Lift was specifically designed to clear walls 10 inches thick and up to 58 inches high with push button ease, making it perfect for your above ground backyard pool or spa.

All stainless steel construction with the same durable powder coating provided on all of our other lifts.

Features:

- 400 Pound (181 kg) weight capacity
- Rechargeable battery operation with charger
- Clears walls up to 58" high, 10" thick
- Removable manual rotation bar
- Durable mesh sling- standard size: 38"x20"
- Adjustable stainless steel chains
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty

Options:

- *Hard seat option - Item #F-0320HSO*
- *Youth sling - Item #F-34EZCS / 34"x18" (86x46cm)*
- *Extra long sling - Item #F-037EZLS / 50"x28" (127x71cm)*
- *Standard sling with rigid insert & seat belt - Item #F-036RSS*
- *Additional standard in ground anchor - Item #F-612TLA*
- *Cover - Item #F-440HBC*
- *Upgrade Pack - Item #F-03PEZAP - Includes battery, rigid sling with seat belt, sling bars and cover*



Sling with rigid seat and belt



Hard Seat Option

PROUDLY SOLD BY:



104 E. CALTON RD. SUITE 110
 LAREDO, TX 78041
 PH. (956) 723-0770

WWW.SPECIALMOBILITY.COM