



Aqua Creek Products

The Total Solution for Pool and Spa Access



The Titan 600™ Lift

Item# F-600TL

**600 lb.
capacity**

The Titan 600™ Lift shares a lot of the same great features as the Revolution Lift but is designed to accommodate patrons up to 600 lbs. No other ADA compliant lift on the market has this weight capacity. **Now independently verified by Wyle Labs to be ADA compliant.**

This lift provides the ability to service multiple bodies of water as well as both above ground spas and in ground pools and spas.

Features:

- 600 Pound (272 kg) weight capacity
- Independently verified to meet ADA standards
- Rechargeable battery with wall mount charger
- Clears 30" high walls, up to 10" thick
- 360 degree powered rotation in either direction
- Adjustable footrest & lap belt
- Dual flip-up arm rests
- Simple push button operation
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty

Options:

- *Upgrade Pack - Item #F-600TLAP*
Includes battery, headrest, chest strap & cover
- *Headrest - Item #F-722HRR*
- *Chest strap - Item #F-423CS*
- *Protective Cover - Item #F-720RLC*
- *Transport Cart - Item #F-614TTC*
- *Cycle Attachment - Item #F-019CA*



Inspected And Verified to Meet
Applicable Sections of ADA
Guidelines Section 1009.2 By:

wyle
laboratories

A Nationally Recognized
Test Laboratory



PROUDLY SOLD BY:



104 E. CALTON RD. SUITE 110
LAREDO, TX 78041
PH. (956) 723-0770
WWW.SPECIALMOBILITY.COM