

# CHARIOT™

Patent Pending



***It's a Lift on Wheels!***



*Transport your scooter or  
powerchair safely with a small  
or mid-size car. Incredible!*

 **BRUNO®**  
for your independence

# It's a Lift on Wheels!



## A Most Remarkable Development!

Many people who depend on scooters and power wheelchairs for mobility have reluctantly traded in their small or mid-size cars for vans, minivans or SUVs to transport their personal mobility devices.

Now there's an alternative! Introducing the Chariot – an incredible new “lift on wheels” that empowers a small vehicle to handle a load of up to 350 lb with ease.

Positioned behind a car equipped with a frame mounted Class I or II receiver hitch, the Chariot uses two independent wheels that swivel – following the car precisely. With its own suspension, this remarkable lift virtually eliminates any additional weight or force on the car!

It's a miracle!

## Look At What People Are Saying About the Chariot...


*“It's great! I don't feel the vehicle bounce like with other common lifts.”*

*“It's so easy to drive...the lift follows right along, even in reverse.”*


*“I think it's very nice. You can't even feel it back there.”*

*“It's not like a trailer..so simple to maneuver!”*


*“It's just like an extension of my car!”*




# The Chariot Works On Many Different Cars!




When used on a small car, it's shorter than a Full-Size Pickup Truck



Minivan



Chariot On Small Car



Full-Size Pickup Truck




**Keep All Your Vehicle Seats!**



**Fits In Most Garages!**

**The Chariot Goes Where You Go... Just Drive Normally!**



**NO JACKKNIFING**

*Swivel wheels on the Chariot means the lift follows you precisely – No “Jackknifing” – You might not even know it’s there.*

**The Chariot Does Not Change Position!**

# Chariot Features...



Self-Tensioning Retractable Securement Belts - One Hand Operation



Swivel Wheels With Independent Suspension



Easy Fold Platform



Drive On/Off Platform From Either Side - No Reversing



Self-Leveling Platform



Optional Corner Markers



Running/Brake Lights and Turn Signals



Manual Back-Up



Lighted License Plate Holder (for use in some states/provinces)



One Button Operation with On/Off Key



Three or Four Wheeled Scooters



Mid-Wheel Drive Powerchairs

The Chariot Works With Many Types Of Mobility Devices...



Front-Wheel/Rear-Wheel Drive Powerchairs



Manual Wheelchairs with Seating Systems

## KEY FEATURES


	CHARIOT
Maximum Weight Capacity	350 lbs/159 kg
Power Up/Down	yes
Manual Assist Fold	yes
Single Switch Operation	yes
Adjustable Platform	yes
Powered By Vehicle Battery	yes
Emergency Backup	yes
Compatible w/Class I or Class II receiver hitch (1 1/4" square) - frame mount only	yes
Uses Flat-4 Wiring Connector (Required)	yes
Optional Spare Tire Available	yes
3 Year Limited Warranty	yes
Can Reinstall on Another Vehicle*	yes
Eligible for Auto Manufacturer Rebate*	yes

**\*IMPORTANT:** Contact your authorized Bruno dealer for specific details on your vehicle!

## Other great Bruno products...



ELITE STAIRLIFT



CUSTOM CURVED RAIL STAIRLIFT



VERTICAL PLATFORM LIFT



CURB-SIDER® VEHICLE LIFT

**NOTE:** The Chariot was designed and tested for use with sedans using a frame-mounted Class I or II receiver hitch. This product is not intended for use with other vehicles such as minivans, SUVs, pickup trucks, etc.

## See The Chariot In Action...



[www.bruno.com](http://www.bruno.com)



OUT-SIDER® MERIDIAN™ VEHICLE LIFT



VALET™ SIGNATURE SEATING

Represented by:

**BRUNO INDEPENDENT LIVING AIDS, INC.®**

1780 Executive Drive, P.O. Box 84, Oconomowoc, WI 53066

[www.bruno.com](http://www.bruno.com)

