

Valet™ Limited

Enter and exit your vehicle with ease!



Shown on: Toyota RAV4



This marvelous seat rotates, then automatically reclines as it passes through the door frame for **extra head clearance.**



Seat Rotates



Comes Out and Lowers



Shown on: Chrysler Town & Country

BRUNO INDEPENDENT LIVING AIDS, INC.®

www.bruno.com



Valet™ Limited Features



Pre-Programmed for maximum legroom



Easy one button operation



Integrated Flip-up Footrest



Standard Valet Limited black fabric seat



Seat-Mounted Controls

- 1.) Power Recline – Both Sides of Seat
- 2.) Power Forward/Backward Seat Adjustment – Both Sides of Seat
- 3.) Unit Operation – Outboard Side Only

Seat Color and Material Options



See material swatches for accurate color representation.



Katzkin leather option handsomely complements existing factory seats.



Vehicle color-matched Katzkin leather for select Toyota/Lexus models*

*Requires interior color code of the vehicle

APPLICATIONS

| Model | Year Range |
|--------------------------|------------|
| Chrysler Town & Country | 2008-12 |
| Dodge Grand Caravan | 2008-12 |
| Ford Escape/Hybrid | 2001-12 |
| Honda CR-V | 2007-12 |
| Honda Odyssey | 2005-10 |
| Lexus RX330 | 2004-06 |
| Lexus RX350 | 2007-12 |
| Lexus RX400/450h | 2006-12 |
| Lexus LX570 | 2008-12 |
| Mercury Mariner/Hybrid | 2001-11 |
| Mazda Tribute/Hybrid | 2001-12 |
| Toyota Highlander/Hybrid | 2004-12 |
| Toyota RAV4 | 2006-12 |
| Toyota Sienna | 2004-12 |
| Toyota Venza | 2009-12 |

OTHER FEATURES

| | |
|---|----------------|
| Maximum Weight Capacity | 330 lb/ 150 kg |
| Eligible for Auto Manufacturer Rebate | yes |
| Meets All Applicable Federal Motor Vehicle Safety Standards (FMVSS) | yes |
| Uses Automotive Seat Belt | yes |
| 3 Year Warranty | yes |
| Requires Structural Modification to Vehicle | no |

IMPORTANT: Contact your authorized Bruno dealer for specific details on your vehicle!