

Valet™

Signature Seating



Bruno's Revolutionary Seating

Introducing...

Valet™
Signature Seating

BRUNO
for your independence

Valet
Valet Plus
Valet Limited
Valet LV



Valet Signature Seating

Bruno Independent Living Aids, the world's leader in vehicle lifts and accessible seating for the mobility marketplace, has launched a complete line of turning seats — Valet Signature Seating. Acknowledging that everyone deserves a bit of luxury, Valet Signature Seats deliver stylish comfort and elegant assistance.

A wide variety of vehicles can accommodate a Valet Signature Seat, and every seat rewards its user with reliable performance and effortless access into and out of their vehicle.



Simple, one button operation is standard for all Valet Signature Seats, making them extremely easy to use. Nothing has been sacrificed to deliver maximum safety. Every model has been thoroughly tested by independent, licensed agencies. And, in most cases, the use of the original vehicle seat belt is retained for ultimate safety.

Valet Signature Seating will present an entirely new opportunity for your life on the road. Grab it!

Valet

A Workhorse...

A simple and reliable solution for access to higher vehicles, such as today's popular minivans, SUVs, CUVs and pickup trucks, the Valet Seat provides full power rotation and lowering to eliminate the twisting and turning associated with entering and exiting many vehicles.

Many people like to ride above the traffic with a nice unobstructed view. Now, you can enjoy the ride with a seat that gracefully places you in control. The Valet requires no structural vehicle modification, thus allowing reinstallation of the original seat, sustaining the vehicle's value for resale.



Available in a broad range of colors, in two upholstery options - Ultraleather™ Plus and Metro Tech fabric.

Assists in entering and exiting your higher vehicle.

Easy one button operation hand-held pendant

Valet Plus

An Elegant Addition...

A step-up from the Valet, the Valet Plus is a luxury seat that combines style with marvelous functionality, the Valet Plus represents the reward you deserve. Individually programmed at installation to ensure the best fit and function, the Valet Plus features power forward/backward seat adjustment for maximum legroom, and full seat recline. An integrated flip-up footrest and retention of the original seat belt provides both comfort and security.

Designed with factory aesthetics in mind, the Valet Plus is available in an impressive array of colors and two material choices - Metro Tech fabric and Ultraleather™ Plus. It feels like the seat has been tailored just for you!



Full recline with levers on both sides



A smaller footprint and lower unit profile means more legroom for first and second row passengers.



Power forward/backward seat adjustment



For more information on this seat please refer to the Valet Plus brochure!

Valet Limited

A Stylish Advantage...

An engineering marvel, the Valet Limited seat provides an array of powered functions: power forward/backward seat adjustment; power rotation/raising and lowering outside the vehicle; and power seat recline. The seat is even programmed to automatically recline as it passes through the door opening, providing extra headroom. Good looks come standard on this premium seat as it “disappears” into the interior of select vehicles. Black pattern fabric seat material is standard with a variety of leather and fabric options available.

Katzkin®

Vehicle color-matched Katzkin leather for select Toyota/Lexus models



The seat rotates...



...then automatically reclines to provide head clearance...



...as it comes out of the vehicle and lowers to the desired height!



For more information on this seat please refer to the Valet Limited brochure!

Valet LV

A Convenient Seat...

The simplest of all the Valet Signature Seating solutions, the Valet LV is available for both the driver (front left) and passenger (front right) positions. It is ideal for assistance in getting in and out of lower vehicles. At the touch of a button, the Valet LV rotates to provide easy access to the seat. Caregivers marvel that the Valet LV rotates out over the door sill, with a manual “seat slider” that bridges the gap for an easier transfer.

You don't have to drive a minivan or SUV to enjoy the magic of Valet Signature Seating. The Valet LV is perfect for sedans and other lower vehicles. Try it for yourself!



Caregivers marvel at how the Valet LV rotates the seat over the door sill, bridging the gap for an easier transfer.



At the touch of a button, the seat power rotates to provide easy access to the seat.



The Valet LV can handle a user weight up to 330 lbs (159 kg).

Metro Tech Fabric



Ultraleather™ Plus



See material swatches for accurate color representation.

Comfort Seat



A Revolutionary Offering in Automotive Seating

Crafted with the highest polyurethane resins on the market, only Ultraleather™ Plus and Metro Tech fabric give you softness, comfort and luxury while offering unmatched durability to protect against the bumps, spills and scrapes of life on the road. These materials offer many benefits that original vehicle seats do not.

New Comfort Seat Highlights

- OEM aesthetics
- State of the art foam technology
- Easy recline levers on both sides
- Earth friendly (Greenguard® certified)
- Many color choices
- Matching armrests available
- Easy to clean – simply wipe up stains with soap and water
- Disinfects with a 5:1 water/bleach solution
- Micro-foam layer under upholstery for long-lasting comfort
- Controls temperature for skin protection
- Scratch and blemish resistant – pet and kid friendly
- Meets all applicable FMVSS requirements



Attention Consumers! Why it's unwise to use an OEM factory seat:

Using an original vehicle seat in a turning application can create serious problems:

- The rider may be outside of the "intended position" set by the vehicle manufacturer, which may result in seat belts and airbags not functioning properly, as well as unsafe head clearance and improper visibility
- This will add significant installation fees when modifying the original vehicle seat to make it fit on the base unit
- It could be expensive or impossible to return the vehicle to its original condition for resale

Here are just a few other great Bruno products . . .



ELAN STAIRLIFT



VERTICAL PLATFORM LIFT



CUSTOM CURVED RAIL STAIRLIFT



CURB-SIDER®



JOEY™ VEHICLE LIFT



CHARIOT™ VEHICLE LIFT

Represented by:

Innovation . . . for your Independence & Convenience

KEY FEATURES

	Valet	Valet Plus	Valet Limited	Valet LV
Maximum Weight Capacity	330 lb/ 150 kg	370 lb/ 170 kg	330 lbs/ 150 kgs	330 lbs/ 150 kgs
Power Up/Down	yes	yes	yes	N/A
Power Rotate	yes	yes	yes	yes
3 Year Warranty	yes	yes	yes	yes
Hand-Held Pendant Controls	yes	yes	yes	yes
Powered by Vehicle Battery	yes	yes	yes	yes
Requires Structural Modification to Vehicle	no	no	no	no
Can Reinstall in Another Vehicle	yes	yes	yes	yes
Eligible for Auto Manufacturer Rebate	yes	yes	yes	yes
Crash Tested	yes	yes	yes	yes
Uses Automotive Seat Belt	yes*	yes	yes	yes

IMPORTANT: Contact your authorized Bruno dealer for specific details on your vehicle!

*Some applications require the Bruno EXO-Frame seat belt replacement system

BRUNO INDEPENDENT LIVING AIDS, INC.®

1780 Executive Drive, Oconomowoc, WI 53066 • www.bruno.com

