



Jazzy[®]

1450

The Jazzy[®] 1450 is a front-wheel drive bariatric power base which has a weight capacity of 600 lbs. and durable construction to make it an exceptionally rugged power chair. Its unique design delivers high-performance operation indoors and out.



Pride[®]
Mobility Products Corp.
LIVE YOUR BEST[®]

SPECIFICATIONS

Weight capacity	600 lbs.
Maximum speed	Up to 4.5 mph
Ground clearance	3" at center frame
Turning radius	26.5"
Overall length	38.25" without front riggings
Overall width	29.25"
Drive wheels	14" solid
Front wheels	4"
Rear wheels	9"
Suspension	Swing Arm Suspension
Drivetrain	Two-motor, front-wheel
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	90 Amp VR2
Specialty controls	N/A
Battery charger	Off-board, 8 amp
Per-charge range	Up to 12 miles
Battery requirements	Two 12 volt, deep cycle, Group 24, 70 Ah
Battery weight	52 lbs. each
Base weight	195 lbs.
Standard seat weight	37 lbs. (medium-back)
Max. seat dimensions	32" x 28"
Warranty	Lifetime limited warranty on frame; 2-year warranty on electronics; 18-month warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

MODELS

Jazzy 1450-2VHD-SS	(Code: K0826)
Jazzy 1450-2VHD-C	(Code: K0827)
Jazzy 1450-2SPVHD-SS	(Code: K0839)

FEATURES

- 600 lbs. weight capacity
- Maximum speed up to 4.5 mph
- High-performance, high-efficiency motor package
- Seat dimensions up to 32" x 28"

OPTIONS

- Synergy® Seat
- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

COLOR CHOICES



Jazzy Red



Jazzy Blue