



Jazzy[®]

ELITE HD

The Jazzy[®] Elite HD offers heavy-duty construction, a 450 lbs. weight capacity, and excellent performance in a stylish, highly maneuverable package.



Pride[®]
Mobility Products Corp.
LIVE YOUR BEST[®]

SPECIFICATIONS

Weight capacity	450 lbs.*
Maximum speed	4.2 mph
Ground clearance	3"
Turning radius	24"
Overall length	42.5"
Overall width	24.5"
Drive wheels	14" solid
Front wheels	3" solid anti-tips
Rear wheels	6" solid casters
Drivetrain	Two-motor, front-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	60 amp controller
Specialty controls	N/A
Battery charger	Off-board, 5 amp (standard)
Per-charge range	9.8 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	NF-22
Battery weight	37 lbs. each
Base weight	89 lbs.
Standard seat weight	51.5 lbs. (high-back)
Max. seat dimensions	24" x 24"
Warranty	5-year limited warranty on frame; 13-month limited warranty on electronics; 13-month limited warranty on drive motors

* Maximum weight capacity is reduced to 400 lbs. (181 kg) on 18" deep seat configurations.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

MODELS

Jazzy Elite HD 2HD-SS	(Code: K0824)
Jazzy Elite HD-2HD-C	(Code: K0825)

FEATURES

- Two-motor, front-wheel drive
- 14" knobby tires and optimum weight distribution provide superior traction
- NF-22 batteries for maximum range
- 5 amp, off-board charger

OPTIONS

- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

RECOMMENDED LIFTS

Pride Lifts & Ramps

- Backpacker® Plus
- Outlander
- Backpacker® MV
- PrideBoom 250



COLOR CHOICES



Red



Blue

PRIDE[®]
PERFORMANCE MOBILITY
— LIVE YOUR BEST[®] —

(US) 800-800-8586 ▪ (Canada) 888-570-1113 ▪ www.pridemobility.com