



# Jazzy<sup>®</sup>

## SELECT 6

The Jazzy Select<sup>®</sup> 6 provides maximum front stability combined with the advantages of in-line motor technology. Active-Trac<sup>®</sup> Suspension and 6" front and rear casters offer superior performance indoors and out. An impressive selection of standard convenience features make it the easiest power chair to use.



**Pride<sup>®</sup>**  
Mobility Products Corp.

— LIVE YOUR BEST<sup>®</sup> —

## SPECIFICATIONS

<b>Weight capacity</b>	300 lbs.
<b>Maximum speed</b>	Up to 4 mph
<b>Ground clearance</b>	2"
<b>Turning radius</b>	20.5"
<b>Overall length</b>	36.625" without front riggings
<b>Overall width</b>	23.5"
<b>Wheels</b>	Drive: 10"
<b>Front casters</b>	6" solid
<b>Rear casters</b>	6" solid
<b>Suspension</b>	Active-Trac®
<b>Drivetrain</b>	Two-motor, in-line, mid-wheel drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	40 amp, PG GC 2 controller
<b>Specialty controls</b>	N/A
<b>Battery charger</b>	Off-board, 3 amp
<b>Per-charge range</b>	Up to 15 miles
<b>Battery requirements</b>	(2) 12 volt, deep cycle
<b>Battery size</b>	U-1
<b>Battery weight</b>	24.5 lbs. each
<b>Base weight</b>	100 lbs.
<b>Standard seat weight</b>	37 lbs. (medium-back)
<b>Max. seat dimensions</b>	20" x 20"
<b>Warranty</b>	5-year limited warranty on frame; 13-month limited warranty on electronics; 13-month limited warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

## MODELS

Jazzy Select 6-SS | (Code: K0822)  
Jazzy Select 6-C | (Code: K0823)

## FEATURES

- 6" casters on front & rear for maximum stability
- In-line motor technology for increased efficiency, torque and range
- Articulating rear caster beam
- Easy access charger port and circuit breaker
- Height and depth-adjustable foot platform
- Built with ease of service in mind

## OPTIONS

- Power elevating seat
- Depth-adjustable, fold-flat seat provides fast, simple adjustment for enhanced, personalized comfort
- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

## RECOMMENDED LIFTS

### ***Pride™ Lifts & Ramps***

- Backpacker® MV
- Backpacker® Plus
- Outlander
- SilverBoom 250



## COLOR CHOICES



Red



Blue

**PRIDE**<sup>®</sup>  
PERFORMANCE MOBILITY  
— LIVE YOUR BEST<sup>®</sup> —

(US) 800-800-8586 • (Canada) 888-570-1113 • [www.pridemobility.com](http://www.pridemobility.com)