



Jazzy[®]

SELECT ELITE

The Jazzy Select[®] Elite delivers a reliable blend of power, performance and style. In-line, front-wheel drive technology gives the Elite excellent stability and maneuverability for solid performance indoors and out.



Pride[®]
Mobility Products Corp.
— LIVE YOUR BEST[®] —

SPECIFICATIONS

Weight capacity	300 lbs.
Maximum speed	Up to 4 mph
Ground clearance	1.5" at anti-tips
Turning radius	24.75"
Overall length	34.375" without front riggings
Overall width	22.75"
Drive wheels	9"
Front casters	3" solid
Rear casters	6" solid
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG GC 3 controller
Specialty controls	N/A
Battery charger	Off-board, 3 amp
Per-charge range	Up to 15 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	U-1
Battery weight	24.5 lbs. each
Base weight	68.5 lbs.
Standard seat weight	43.5 lbs. (comfort high-back)
Max. seat dimensions	20" x 20"
Warranty	5-year limited warranty on frame; 13-month limited warranty on electronics; 13-month limited warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

MODELS

Jazzy Select Elite-2S-SS	(Code: K0822)
Jazzy Select Elite-2S-C	(Code: K0823)

FEATURES

- In-line, front-wheel drive technology
- 40 amp, PG GC 3 controller
- Easy access to batteries from rear
- Simple main frame design for easy serviceability
- New comfort high-back seat with headrest
- New Jazzy® Armor™ includes shroud and controller guards to protect from daily wear and tear

OPTIONS

- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Deluxe Jazzy® Armor™ package (coming soon)

RECOMMENDED LIFTS

Silver Star®
Lifts & Ramps

- Backpacker® Plus
- Outlander
- Backpacker® MV



COLOR CHOICES



Jazzy Red



Jazzy Blue

PRIDE®
PERFORMANCE MOBILITY
— LIVE YOUR BEST® —

(US) 800-800-8586 • (Canada) 888-570-1113 • www.pridemobility.com