

THE ULTIMATE
HIGH-PERFORMANCE
10" SCOOTER

VICTORY® SPORT

The stylish, high-performing Victory Sport Scooter features innovative design and a full complement of unique standard features. A wrap-around delta tiller, high-visibility LED lighting and full suspension are just a few features that make the Victory Sport an exceptional scooter choice.



Victory[®] Sport 4-wheel

Model number:	SC710 DXW
Weight capacity:	350 lbs.
Maximum speed ¹ :	Up to 8 mph
Ground clearance ² :	2.25"
Turning radius:	51.75"
Overall length ² :	47"
Overall width ² :	26"
Range per charge ^{1,5} :	Up to 11 miles
Tire size:	Front: 3.5x10 solid Rear: 4x10 solid
Total weight without batteries ^{1,2} :	150 lbs.
Heaviest piece when disassembled ^{2,6} :	Front: 68 lbs.
Standard seating:	Type: Deluxe reclining high-back with sliders and headrest Weight: 54 lbs. Material: Black Vinyl Dimensions: 18" width 18" depth
Seat-to-deck range ² :	17.25" - 18.75"
Drivetrain:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual braking system:	Electronic, regenerative and electromechanical
Battery requirements ³ :	Two 12-volt, deep cycle Size: 40 AH Weight ⁴ : 32.5 lbs. each
Battery charger:	5-amp, off-board
Warranty:	3-year limited on frame 3-year pro-rata on drivetrain 1-year limited on electronics

Colors:



Candy Apple Red



Viper Blue

¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

²Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³AGM or gel cell type required.

⁴Battery weight may vary based on manufacturer.

⁵Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.

⁶Includes standard seat weight.

NOTE: All specifications subject to change without notice.

Features

- Exclusive one-hand feather-touch disassembly provides the ultimate quick transport system
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- Exclusive stylish, lightweight, non-scuffing, black, low-profile wheels
- Backlit battery gauge
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Lightweight seat with viscoelastic foam inserts for increased comfort and sliders
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port

Options

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag

Pride[™] Lifts & Ramps

- Backpacker[®] Plus
- Outlander
- Backpacker[®] MV



Backpacker[®] Plus



Outlander

These lifts are compatible with the Victory[®] Sport.