

SHOPRIDER



K0823

6RUNNER[®] 10 [888WNLM]

Much like our original 6Runner 10 power chair, the NEW generation 6Runner 10 incorporates enhanced stability having all 6 points of contact with the ground, and true mid-wheel technology for maneuverability that matches a more natural feeling for an independent lifestyle. Choosing the NEW 6Runner 10 means choosing a luxury product at an economy price!

Toll-Free: 800.743.0772
www.shoprider.com



swing away option joystick
on 6Runner 10 with VR2 control

MID WHEEL DRIVE

- 300lb weight capacity
- Height adjustable seat & adjustable arms
- Stable Four post seating system
- Full suspension
- 18" wide high back recline executive seat
- Easy Battery Access
- 18" wide seat and other seat options available

6RUNNER® 10



Overall Dimensions (L x W x H)	in	42 x 23 x 46
Number, Size of Tires (front; rear)	Front (rear)	2, 7"; 2, 10"; 2, 7"
Suggested User Weight (on a level road)	lb	300
Total Weight (w/35 Ah batteries, seat)	lb	191
Each Battery Module Weight	lb	5.8
Battery Module Capacity	_V_Ah x Pcs	12V35Ah x 2
Charger		Off Board (5A)
Output Power of Motor	hp	0.6 x2
Maximum Speed ²	mph	5
Max. Climbing Angle/Safe Climbing Angle(based on user weight of 90kg/200lb.)	% (°)	10 (6)
Range (per charge) (After the battery and mechanical moving parts fully break in)	mile	16.2
Turning Radius	in	15.75

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant

b) Ground surface

c) Battery capacity and conditions

d) Type of charger

e) Ambient temperature

f) The way of driving

g) If the battery and mechanical moving parts fully break in

h) Etc.