


SHOPRIDER



K0822

FPC [PHFW-III8, PHFW-II20]

Presenting Shoprider's folding power chair, the FPC. The FPC is packed with pioneering technology and style. The steel black frame, puncture proof tires and height adjustable arms are all standard with the FPC. The standard seat cushion on the FPC makes it a pleasure to drive.



Seat Cushion Option

Toll-Free: 800.743.0772
www.shoprider.com



Adjustable Height Arm Rest



Flip-Up Arm Rest



Adjustable Length VSI Controller

FOLDING POWERCHAIR

- Optional Seat and back cushion for supreme comfort
- Available in 18" and 20" seat base
- Optional elevating leg rests
- Adjustable length VSI controller
- Easy to fold for transportation

FPC



Black

		PHFW-III8	PHFW-II20
Overall Dimensions (L x W x H)	in	46 x 25 x 37	46 x 27 x 37
Number, Size of Tires (front; rear)	Front (rear)	2, 8"; 2, 13"	2, 8"; 2, 13"
Suggested User Weight (on a level road)	lb	300	300
Total Weight (w/U1*2)	lb	80	80
Each Battery Module Weight (U1)	lb	25 ea x2	25 ea x2
Battery Module Capacity	_V_Ah x Pcs	12V36Ah x 2	12V36Ah x 2
Charger		Off Board (4A)	Off Board (4A)
Output Power of Motor	hp	0.6	0.6
Maximum Speed ²	mph	5	5
Max. Climbing Angle/Safe Climbing Angle(based on user weight of 90kg/200lb.)	% (°)	10 (6)	10 (6)
Ground Clearance	in	3	3
Turning Radius	in	24	24

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant

b) Ground surface

c) Battery capacity and conditions

d) Type of charger

e) Ambient temperature

f) The way of driving

g) If the battery and mechanical moving parts fully break in

h) Etc.